

Product Datasheet

Human Vascular Endothelial Growth Factor D (VEGFD) ELISA Kit (orb1209018)

Description	Human Vascular Endothelial Growth Factor D (VEGFD) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O43915
Sensitivity	6.1pg/mL
Entrez	2277
Expiration Date	Please enquire.