

## Product Datasheet

# Human Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit (orb1209037)

<b>Description</b>	Human Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O14625</a>
<b>Sensitivity</b>	10.8pg/mL
<b>Entrez</b>	<a href="#">6373</a>
<b>Expiration Date</b>	Please enquire.