

Product Datasheet

Canine Interleukin 8 (IL8) ELISA Kit (orb1209043)

| | |
|--------------------------|--|
| Description | Canine Interleukin 8 (IL8) ELISA Kit |
| Reactivity | Canine |
| Range | 0.156-10ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | P41324 |
| Sensitivity | 0.058ng/mL |
| Entrez | 403850 |
| Expiration Date | Please enquire. |