

## Product Datasheet

### Guinea pig Albumin (ALB) ELISA Kit (orb1209053)

<b>Description</b>	Guinea pig Albumin (ALB) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	1.56-100µg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q6WDN9</a>
<b>Sensitivity</b>	0.48µg/mL
<b>Entrez</b>	<a href="#">100379232</a>
<b>Expiration Date</b>	Please enquire.