

## Product Datasheet

# Human Enterovirus 71 Antibody IgG ELISA Kit (orb1209091)

<b>Description</b>	Human Enterovirus 71 Antibody IgG, EV71 AB IgG ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	EV71 ab IgG
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Qualitative
<b>Sample Types</b>	Serum, plasma, cell culture supernatant
<b>Expiration Date</b>	Please enquire.