

Product Datasheet

Rat Hydroxyproline ELISA Kit (orb1209102)

Description	Rat Hydroxyproline, HYP ELISA Kit
Reactivity	Rat
Range	0.05 - 20ng/ml
Tested Applications	ELISA
Target	Hyp
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	25.6ng/mL
Sensitivity	0.027ng/ml
Expiration Date	Please enquire.