

Product Datasheet

Rat Adrenocorticotrophic Hormone ELISA Kit (orb1209109)

Description	Rat Adrenocorticotrophic Hormone, ACTH ELISA Kit
Reactivity	Rat
Range	0.3 - 90ng/L
Tested Applications	ELISA
Target	ACTH
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	96ng/L
Sensitivity	0.13ng/L
Expiration Date	Please enquire.