

Product Datasheet

Human 4-hydroxynonenal ELISA Kit (orb1209121)

Description	Human 4-hydroxynonenal, 4HNE ELISA Kit
Reactivity	Human
Range	10 - 2000ng/L
Tested Applications	ELISA
Target	4HNE
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	2400ng/L
Sensitivity	5.23ng/L
Expiration Date	Please enquire.