

## Product Datasheet

### Fish Estradiol ELISA Kit (orb1209131)

<b>Description</b>	Fish Estradiol, E2 ELISA Kit
<b>Reactivity</b>	Fish
<b>Range</b>	15 - 3000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	E2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive
<b>Sample Types</b>	Serum, plasma, cell culture supernatant
<b>Solubility (25°C)</b>	3200ng/L
<b>Sensitivity</b>	7.51ng/L
<b>Expiration Date</b>	Please enquire.