

Product Datasheet

Human Calsyntenin 1 ELISA Kit (orb1209153)

Description	Human Calsyntenin 1, CLSTN1 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CLSTN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	109,793
Uniprot ID	O94985
Solubility (25°C)	4000ng/L
Sensitivity	12.32ng/L
Entrez	22883
Expiration Date	Please enquire.