

## Product Datasheet

# Human Platelet-derived Growth Factor Receptor Like Protein ELISA Kit (orb1209164)

<b>Description</b>	Human Platelet-derived Growth Factor Receptor Like Protein, PDGFRL ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-30ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	PDGFRL
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	41,861
<b>Uniprot ID</b>	<a href="#">Q15198</a>
<b>Solubility (25°C)</b>	32ng/ml
<b>Sensitivity</b>	0.028ng/ml
<b>Entrez</b>	<a href="#">5157</a>
<b>Expiration Date</b>	Please enquire.