

Product Datasheet

Human Nadh-ubiquinone Oxidoreductase 75 kDa Subunit ELISA Kit (orb1209196)

Description	Human Nadh-ubiquinone Oxidoreductase 75 kDa Subunit, NDUFS1 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	NDUFS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	79,468
Uniprot ID	P28331
Solubility (25°C)	4000ng/L
Sensitivity	11.11ng/L
Entrez	4719
Expiration Date	Please enquire.