

Product Datasheet

Human Peptidyl-prolyl Cis-Trans Isomerase D ELISA Kit (orb1209212)

Description	Human Peptidyl-prolyl Cis-Trans Isomerase D, CYCLOPHILIN-40 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	Cyclophilin-40
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,764
Uniprot ID	Q08752
Solubility (25°C)	16ng/ml
Sensitivity	0.013ng/ml
Entrez	5481
Expiration Date	Please enquire.