

Product Datasheet

Human Bdnf/Nt-3 Growth Factors Receptor ELISA Kit (orb1209230)

Description	Human Bdnf/Nt-3 Growth Factors Receptor, NTRK2 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	NTRK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	91,999
Uniprot ID	Q16620
Solubility (25°C)	96ng/ml
Sensitivity	0.13ng/ml
Entrez	4915
Expiration Date	Please enquire.