

Product Datasheet

Human Regulatory-associated Protein of Mtor ELISA Kit (orb1209235)

Description	Human Regulatory-associated Protein of Mtor, RPTOR ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	Raptor
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	149,038
Uniprot ID	Q8N122
Solubility (25°C)	2400ng/L
Sensitivity	5.18ng/L
Entrez	57521
Expiration Date	Please enquire.