

Product Datasheet

Human C5A Anaphylatoxin Chemotactic Receptor ELISA Kit (orb1209259)

Description	Human C5A Anaphylatoxin Chemotactic Receptor, C5AR1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	C5AR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,864
Uniprot ID	Q16581
Solubility (25°C)	4800ng/L
Sensitivity	9.77ng/L
Entrez	719
Expiration Date	Please enquire.