

## Product Datasheet

# Human B- and T-lymphocyte Attenuator ELISA Kit (orb1209285)

<b>Description</b>	Human B- and T-lymphocyte Attenuator, BTLA ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.2-60ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	BTLA
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	32,834
<b>Uniprot ID</b>	<a href="#">Q7Z6A9</a>
<b>Solubility (25°C)</b>	64ng/ml
<b>Sensitivity</b>	0.095ng/ml
<b>Entrez</b>	<a href="#">151888</a>
<b>Expiration Date</b>	Please enquire.