

Product Datasheet

Human Vitamin K-dependent Protein Z ELISA Kit (orb1209323)

Description	Human Vitamin K-dependent Protein Z, PROZ ELISA Kit
Reactivity	Human
Range	70-15000ng/mL
Tested Applications	ELISA
Target	PROZ
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,744
Uniprot ID	P22891
Solubility (25°C)	16000ng/mL
Sensitivity	34.39ng/mL
Entrez	8858
Expiration Date	Please enquire.