

Product Datasheet

Human Uv Excision Repair Protein Rad23 Homolog B ELISA Kit (orb1209329)

Description	Human Uv Excision Repair Protein Rad23 Homolog B, RAD23B ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	RAD23B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,171
Uniprot ID	P54727
Solubility (25°C)	3200ng/L
Sensitivity	7.65ng/L
Entrez	5887
Expiration Date	Please enquire.