

Product Datasheet

Human Ubiquitin-associated and Sh3 Domain-containing Protein A ELISA Kit (orb1209337)

Description	Human Ubiquitin-associated and Sh3 Domain-containing Protein A, UBASH3A ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	UBASH3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	74,123
Uniprot ID	P57075
Solubility (25°C)	2400ng/L
Sensitivity	5.14ng/L
Entrez	53347
Expiration Date	Please enquire.