

Product Datasheet

Human Ubiquitin Thioesterase Otub2 ELISA Kit (orb1209338)

Description	Human Ubiquitin Thioesterase Otub2, OTUB2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	OTUB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,213
Uniprot ID	Q96DC9
Solubility (25°C)	3200ng/L
Sensitivity	7.31ng/L
Entrez	78990
Expiration Date	Please enquire.