

Product Datasheet

Human Tyrosine-protein Phosphatase Non-receptor Type 6 ELISA Kit (orb1209343)

Description	Human Tyrosine-protein Phosphatase Non-receptor Type 6, PTPN6 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	PTPN6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,561
Uniprot ID	P29350
Solubility (25°C)	2400ng/L
Sensitivity	5.18ng/L
Entrez	5777
Expiration Date	Please enquire.