

Product Datasheet

Human SVEP1 ELISA Kit (orb1209402)

Description	Human SVEP1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	SVEP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	390,170
Uniprot ID	Q4LDE5
Solubility (25°C)	3200ng/L
Sensitivity	7.28ng/L
Entrez	79987
Expiration Date	Please enquire.