

## Product Datasheet

# Human Sparc-related Modular Calcium-binding Protein 2 ELISA Kit (orb1209421)

<b>Description</b>	Human Sparc-related Modular Calcium-binding Protein 2, SMOC2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	SMOC2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	49,674
<b>Uniprot ID</b>	<a href="#">Q9H3U7</a>
<b>Solubility (25°C)</b>	24ng/mL
<b>Sensitivity</b>	0.024ng/mL
<b>Entrez</b>	<a href="#">64094</a>
<b>Expiration Date</b>	Please enquire.