

Product Datasheet

Human Sparc-like Protein 1 ELISA Kit (orb1209422)

Description	Human Sparc-like Protein 1, SPARCL1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	SPARCL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	75,208
Uniprot ID	Q14515
Solubility (25°C)	48ng/mL
Sensitivity	0.049ng/mL
Entrez	8404
Expiration Date	Please enquire.