

Product Datasheet

Human Radical S-Adenosyl Methionine Domain-containing Protein 2 ELISA Kit (orb1209489)

Description	Human Radical S-Adenosyl Methionine Domain-containing Protein 2, RSAD2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	RSAD2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,170
Uniprot ID	Q8WXG1
Solubility (25°C)	3200ng/L
Sensitivity	7.65ng/L
Entrez	91543
Expiration Date	Please enquire.