

Product Datasheet

Human Radiation-inducible Immediate-Early Gene Iex-1 ELISA Kit (orb1209490)

Description	Human Radiation-inducible Immediate-Early Gene Iex-1, IER3 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	IER3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,903
Uniprot ID	P46695
Solubility (25°C)	4800ng/L
Sensitivity	9.18ng/L
Entrez	8870
Expiration Date	Please enquire.