

Product Datasheet

Human Protein Wnt-5B ELISA Kit (orb1209510)

Description	Human Protein Wnt-5B, WNT5B ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	WNT5B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,323
Uniprot ID	Q9H1J7
Solubility (25°C)	24ng/mL
Sensitivity	0.03ng/mL
Entrez	81029
Expiration Date	Please enquire.