

Product Datasheet

Human Protein Fam3C ELISA Kit (orb1209524)

Description	Human Protein Fam3C, FAM3C ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	FAM3C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,680
Uniprot ID	Q92520
Solubility (25°C)	9600ng/L
Sensitivity	14.15ng/L
Entrez	10447
Expiration Date	Please enquire.