

Product Datasheet

Human Protein Disulfide-isomerase A4 ELISA Kit (orb1209528)

Description	Human Protein Disulfide-isomerase A4, PDIA4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	PDIA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	72,932
Uniprot ID	P13667
Solubility (25°C)	3200ng/L
Sensitivity	7.58ng/L
Entrez	9601
Expiration Date	Please enquire.