

Product Datasheet

Human Prenylcysteineoxidase 1 ELISA Kit (orb1209545)

Description	Human Prenylcysteineoxidase 1, PCYOX1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	PCYOX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,640
Uniprot ID	Q9UHG3
Solubility (25°C)	4800ng/L
Sensitivity	9.51ng/L
Entrez	51449
Expiration Date	Please enquire.