

Product Datasheet

Human Phosphorylated Microtubule-associated Protein Tau ELISA Kit (orb1209561)

Description	Human Phosphorylated Microtubule-associated Protein Tau, PMAPT ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	pMAPT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	78,928
Uniprot ID	P10636
Solubility (25°C)	4800ng/L
Sensitivity	9.54ng/L
Entrez	4137
Expiration Date	Please enquire.