

Product Datasheet

Human Otoraplin ELISA Kit (orb1209581)

Description	Human Otoraplin, OTOR ELISA Kit
Reactivity	Human
Range	3-900ng/mL
Tested Applications	ELISA
Target	OTOR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,332
Uniprot ID	Q9NRC9
Solubility (25°C)	960ng/mL
Sensitivity	1.64ng/mL
Entrez	56914
Expiration Date	Please enquire.