

Product Datasheet

Human Olfactomedin-like Protein 3 ELISA Kit (orb1209587)

Description	Human Olfactomedin-like Protein 3, OLFML3 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	OLFML3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,010
Uniprot ID	Q9NRN5
Solubility (25°C)	48ng/mL
Sensitivity	0.058ng/mL
Entrez	56944
Expiration Date	Please enquire.