

Product Datasheet

Human Nuclear Receptor Ror-alpha ELISA Kit (orb1209594)

Description	Human Nuclear Receptor Ror-alpha, RORA ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	RORA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,975
Uniprot ID	P35398
Solubility (25°C)	2400ng/L
Sensitivity	4.16ng/L
Entrez	6095
Expiration Date	Please enquire.