

## Product Datasheet

# Human Leucine-rich Repeat Transmembrane Protein Flrt3 ELISA Kit (orb1209681)

<b>Description</b>	Human Leucine-rich Repeat Transmembrane Protein Flrt3, FLRT3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-6000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	FLRT3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	73,004
<b>Uniprot ID</b>	<a href="#">Q9NZU0</a>
<b>Solubility (25°C)</b>	6400ng/L
<b>Sensitivity</b>	9.72ng/L
<b>Entrez</b>	<a href="#">23767</a>
<b>Expiration Date</b>	Please enquire.