

Product Datasheet

Human Interferon-induced Protein with Tetratricopeptide Repeats 1 ELISA Kit (orb1209702)

Description	Human Interferon-induced Protein with Tetratricopeptide Repeats 1, IFIT1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	IFIT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,360
Uniprot ID	P09914
Solubility (25°C)	2560ng/L
Sensitivity	5.16ng/L
Entrez	3434
Expiration Date	Please enquire.