

Product Datasheet

Human Fibronectin Type Iii Domain-containing Protein 3B ELISA Kit (orb1209764)

Description	Human Fibronectin Type Iii Domain-containing Protein 3B, FNDC3B ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	FNDC3B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	132,888
Uniprot ID	Q53EPO
Solubility (25°C)	4800ng/L
Sensitivity	9.59ng/L
Entrez	64778
Expiration Date	Please enquire.