

Product Datasheet

Human Exosome Complex Component Rrp45 ELISA Kit (orb1209770)

Description	Human Exosome Complex Component Rrp45, EXOSC9 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	EXOSC9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,949
Uniprot ID	Q06265
Solubility (25°C)	3200ng/L
Sensitivity	7.91ng/L
Entrez	5393
Expiration Date	Please enquire.