

Product Datasheet

Human Ephrin Type-B Receptor 3 ELISA Kit (orb1209776)

Description	Human Ephrin Type-B Receptor 3, EPHB3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	EPHB3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	110,330
Uniprot ID	P54753
Solubility (25°C)	24ng/ml
Sensitivity	0.027ng/ml
Entrez	2049
Expiration Date	Please enquire.