

Product Datasheet

Human Ephrin Type-B Receptor 2 ELISA Kit (orb1209777)

Description	Human Ephrin Type-B Receptor 2, EPHB2 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	EPHB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	117,493
Uniprot ID	P29323
Solubility (25°C)	96ng/mL
Sensitivity	0.16ng/mL
Entrez	2048
Expiration Date	Please enquire.