

## Product Datasheet

# Human Ephrin Type-B Receptor 1 ELISA Kit (orb1209778)

<b>Description</b>	Human Ephrin Type-B Receptor 1, EPHB1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	30-9000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	EPHB1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	109,885
<b>Uniprot ID</b>	<b>P54762</b>
<b>Solubility (25°C)</b>	9600ng/L
<b>Sensitivity</b>	15.08ng/L
<b>Entrez</b>	<b>2047</b>
<b>Expiration Date</b>	Please enquire.