

Product Datasheet

Human Transmembrane Glycoprotein Nmb ELISA Kit (orb1209834)

Description	Human Transmembrane Glycoprotein Nmb, GPNMB ELISA Kit
Reactivity	Human
Range	0.2-50ng/ml
Tested Applications	ELISA
Target	GPNMB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,923
Uniprot ID	Q14956
Solubility (25°C)	64ng/ml
Sensitivity	0.12ng/ml
Entrez	10457
Expiration Date	Please enquire.