

Product Datasheet

Human Cd177 Antigen ELISA Kit (orb1209851)

Description	Human Cd177 Antigen, CD177 ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	CD177
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,363
Uniprot ID	Q8N6Q3
Solubility (25°C)	128ng/ml
Sensitivity	0.26ng/ml
Entrez	57126
Expiration Date	Please enquire.