

Product Datasheet

Human SLC25A5 ELISA Kit (orb1209875)

Description	Human SLC25A5 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SLC25A5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	32,852
Uniprot ID	P05141
Solubility (25°C)	2400ng/L
Sensitivity	5.47ng/L
Entrez	292
Expiration Date	Please enquire.