

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily ELISA Kit (orb1209932)

Description	Human Tumor Necrosis Factor Receptor Superfamily, Member 9, TNFRSF9 ELISA Kit
Reactivity	Human
Range	2-600ng/L
Tested Applications	ELISA
Target	TNFRSF9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,899
Uniprot ID	Q07011
Solubility (25°C)	800ng/L
Sensitivity	0.96ng/L
Entrez	3604
Expiration Date	Please enquire.