

Product Datasheet

Human Spermine Oxidase ELISA Kit (orb1209941)

| | |
|----------------------------|---|
| Description | Human Spermine Oxidase, SMOX ELISA Kit |
| Reactivity | Human |
| Range | 0.05-30ng/ml |
| Tested Applications | ELISA |
| Target | SMOX |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 61,819 |
| Uniprot ID | Q9NWM0 |
| Solubility (25°C) | 32ng/ml |
| Sensitivity | 0.023ng/ml |
| Entrez | 54498 |
| Expiration Date | Please enquire. |