

Product Datasheet

Human Nectin-4 ELISA Kit (orb1209968)

Description Human Nectin-4, NECTIN4
ELISA Kit

Reactivity Human

Range 10-3800ng/L

Tested Applications ELISA

Target NECTIN4

Storage 2-8°C

Note For research use only

Application notes Sandwich

Sample Types Serum, plasma, cell culture supernates

MW 55,454

Uniprot ID [Q96NY8](#)

Solubility (25°C) 4000ng/L

Sensitivity 5.23ng/L

Entrez [81607](#)

Expiration Date Please enquire.