

Product Datasheet

Human Alcohol Dehydrogenase 6 ELISA Kit (orb1209986)

Description	Human Alcohol Dehydrogenase 6, ADH6 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	ADH6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,073
Uniprot ID	P28332
Solubility (25°C)	24ng/ml
Sensitivity	0.022ng/ml
Entrez	130
Expiration Date	Please enquire.