

## Product Datasheet

# Human Phosphatidylinositol-3,4 ELISA Kit (orb1210032)

<b>Description</b>	Human Phosphatidylinositol-3,4, 5- Trisphosphate 3-phosphatase and Dual-specificity Protein Phosphata...
<b>Reactivity</b>	Human
<b>Range</b>	0.1-40ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	PTEN
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	47,166
<b>Uniprot ID</b>	<b>P60484</b>
<b>Solubility (25°C)</b>	48ng/ml
<b>Sensitivity</b>	0.04ng/ml
<b>Entrez</b>	<b>5728</b>
<b>Expiration Date</b>	Please enquire.